# Service Manual Leica CM3050 S Cryostat

*Leica*

Revision Record

| VERSION<br>DATE | VERSION<br>NUMBER | NEW PAGES     | AMENDMENT DETAIL |
|-----------------|-------------------|---------------|------------------|
| December 2000   | 1.00              | First Version |                  |

# **Service Documentation**

## **Leica CM3050 S**

### **Cryostat**

### **Version 1.0**

Produced by:

Leica Microsystems Nussloch GmbH

Postfach 1120

Heidelberger Str. 17-19

D-69222 Nussloch

Phone : (06224) 143-0

Facsimile : (06224) 143-199

© Leica Microsystems Nussloch GmbH

---

---

## Table of contents

|                  |                       |
|------------------|-----------------------|
| <b>Chapter 1</b> | <b>Introduction</b>   |
| <b>Chapter 2</b> | <b>Mechanics</b>      |
| <b>Chapter 3</b> | <b>Electronic</b>     |
| <b>Chapter 4</b> | <b>Cooling System</b> |
| <b>Chapter 5</b> | <b>Spare Parts</b>    |
| <b>Chapter 6</b> | <b>Other</b>          |

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## Table of contents

|        |                                                                     |    |
|--------|---------------------------------------------------------------------|----|
| 1.1    | General Information .....                                           | 4  |
| 1.1.1  | Dokumentation .....                                                 | 4  |
| 1.1.2  | Specified use and application ....                                  | 5  |
| 1.1.3  | Ordering Spare Parts .....                                          | 6  |
| 1.1.4  | Warranty and Service .....                                          | 7  |
| 1.2    | Safety instructions .....                                           | 8  |
| 1.2.1  | Danger .....                                                        | 8  |
| 1.2.2  | Qualified Personnel .....                                           | 9  |
| 1.2.3  | Symbols used in the text<br>and their meanings .....                | 9  |
| 1.2.4  | Liability .....                                                     | 10 |
| 1.2.5  | Instructions .....                                                  | 11 |
| 1.2.6  | Instrument type .....                                               | 11 |
| 1.2.7  | Transport and installation .....                                    | 11 |
| 1.3    | Instrument-specific information .....                               | 12 |
| 1.3.1  | Information on instrument<br>design and safe handling .....         | 12 |
| 1.3.2  | Integrated safety devices .....                                     | 12 |
| 1.3.3  | Locking the handwheel .....                                         | 13 |
| 1.3.4  | Centering the handwheel grip ..                                     | 14 |
| 1.3.5  | Emergency stop function .....                                       | 14 |
| 1.3.6  | Knife guard .....                                                   | 15 |
| 1.3.7  | Transport .....                                                     | 15 |
| 1.3.8  | Site requirements .....                                             | 15 |
| 1.3.9  | Electrical connections .....                                        | 16 |
| 1.3.10 | Handling microtome<br>knives/blades. ....                           | 16 |
| 1.3.11 | Knife guard/handwheel lock ...                                      | 16 |
| 1.3.12 | Motorized sectioning .....                                          | 17 |
| 1.3.13 | Defrosting/Handling frozen<br>tissue .....                          | 17 |
| 1.3.14 | Frozen parts of the instrument ....<br>and frozen accessories ..... | 17 |
| 1.3.15 | Infectious/<br>radioactive material .....                           | 18 |
| 1.3.16 | Disinfection and cleaning .....                                     | 18 |
| 1.3.17 | Removing/reinstalling the .....                                     | 19 |
| 1.3.18 | Display message<br>'Dry microtome' .....                            | 19 |

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## 1.1 General Information

### 1.1.1 Dokumentation

The information, numerical data, notes and value judgments contained in this manual represent the current state of scientific knowledge and state-of-the-art technology as we understand it following thorough investigation in this field. We are under no obligation to update the present manual periodically and on an ongoing basis according to the latest technical developments, nor to provide our customers with additional copies, updates etc. of this manual.

For erroneous statements, drawings, technical illustrations etc. contained in this manual we exclude liability as far as permissible according to the national legal system applicable in each individual case. In particular, no liability whatsoever is accepted for any financial loss or consequential damage caused by or related to compliance with statements or other information in this manual.

Statements, drawings, illustrations and other information as regards contents or technical details of the present manual are not to be considered as warranted characteristics of our products. These are determined only by the contract provisions agreed between ourselves and our customers.

Leica reserves the right to change technical specifications as well as manufacturing processes without prior notice. Only in this way is it possible to continuously improve the technology and manufacturing techniques used in our products.

This document is protected under copyright laws. Any copyrights of this document are retained by Leica Microsystems Nussloch GmbH.

Any reproduction of text and illustrations (or of any parts thereof) by means of print, photocopy, microfiche, web cam or other methods – including any electronic systems and media – requires express prior permission in writing by Leica Microsystems Nussloch GmbH.

For the instrument serial number and year of manufacture, please refer to the name plate at the back of the instrument.

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© Leica Microsystems Nussloch GmbH

### 1.1.2 Specified use and application

The Leica CM3050 S is a powerful cryostat for routine as well as research applications in biology, medicine and industry.

The instrument has been designed for rapid freezing and sectioning of tissue samples.

The instrument has not been designed for unattended storage of tissue material.

The instrument may only be operated within the scope of its designated use as described above and as per the instructions given in this manual.

Any other use of the instrument is considered improper.

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### 1.1.3 Ordering Spare Parts

To order replacement parts or modules, specify the following information for each part ordered:

1. Product model and serial number.
2. Leica Part number.
3. Part description.
4. Quantity required.

**Spare parts can be obtained from your local Leica office or from:**

Leica Microsystems Nussloch GmbH  
Postfach 1120  
Heidelberger Str. 17-19  
D-69222 Nussloch  
Germany

#### Technical Service - Hotline:

There are two service hotlines for customer service and requests.

Information can be given regarding spare parts, possibilities for repair, service documentation, adjustment and alignment requests etc.

|                      |                                        |
|----------------------|----------------------------------------|
| ☎ : ++ 6224 143 219  | for electronical requests              |
| ☎ : ++ 6224 143 219  | for mechanical requests                |
| Fax: ++ 6224 143 219 | for electrical and mechanical requests |

The agencies are obliged to send complete failure statistics (monthly) to Nussloch for the first year after instruments have been launched onto the market. The results and necessary amendments can be included immediately in the production.

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### 1.1.4 Warranty and Service

Leica Microsystems Nussloch GmbH guarantees that the delivered product has been subjected to a comprehensive quality control procedure based on our strict in-house testing standards in order to ensure that the product complies with its technical specification.

The warranty conditions depend on the contents of the individual contract concluded, supplemented by the warranty conditions of your local Leica sales agency.

Any repairs and/or exchange of parts of the product must be carried out by authorized Leica technical service engineers. Otherwise, any warranty becomes invalid and warranty claims can no longer be made.

The local Leica representative or the manufacturer in Nussloch must be consulted prior to any handling of or changes to the instrument beyond the scope of this instruction manual as well as prior to any modifications or any use of the instrument in combination with non-Leica components not expressly authorized by Leica.

Spare parts and accessories not supplied by Leica can under no circumstances be considered as inspected and/or approved by Leica.

Therefore, installation or use of any such parts may impair the technical design features and thus properties of the instrument.

Leica assumes no liability whatsoever for any damage caused by the use of non-original spare parts or non-original accessories.

The warranty is only valid and warranty claims can only be made as long as the instrument has been operated according to its designated use and according to the instructions given in this manual.

Improper use of the product and/or faulty operation invalidate the warranty and any claims based thereon, and likewise Leica will not assume liability for any consequential damage.

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## 1.2 Safety instructions

### 1.2.1 Danger

While this instrument is in operation, certain components of the unit are inevitably live, carrying a current which can cause severe injuries or death. It is essential to use the precautions mentioned below, in order to reduce the risk of death and/or injury.

1. Only **qualified personnel** (see following page), who are familiar with both the instrument and the instructions supplied together with the instrument, may search the instrument for faults and trouble-shoot and/or repair the instrument.
2. Installation of the instrument must be carried out in compliance with the applicable safety regulations (e.g. DIN, VDE, UL) as well as with any other pertinent national or local rules. Adequate grounding, dimensioning of conductors and corresponding short-circuit protection must be provided, in order to ensure operational safety.
3. During normal operation, all covers have to be installed and must remain closed.
4. Prior to doing visual inspections and/or maintenance work, make sure that the AC power supply is cut off.  
Danger! - prior to cutting off the AC power supply, the instrument is live!
5. If certain measurements have to be done with the current supply on, never touch the electrical connections! The plastic cover at the power supply unit must be installed. Remove all jewelry from your wrists and fingers. Make sure the test devices are in good, operationally safe working order.
6. When working on a live instrument, stand on insulated ground, i.e. make sure that there is no grounding.
7. This list does not necessarily contain all measures that may be necessary for safe operation of the instrument. If you need further information or if specific problems occur, please contact your local Leica office.

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### 1.2.2 Qualified Personnel

as defined in this service manual, are persons who are familiar with installation, maintenance, repair and operation of the product Leica DSC 1, and who have received adequate training in order to carry out their responsibilities. Adequate training is e.g.:

- Professional training, or receiving specific instructions from a professional, or being licensed to connect / disconnect, ground and label electric circuits and instruments / systems according to all applicable safety regulations.
- Professional training, or receiving specific instructions from a professional on how to use / maintain applicable safety regulations, the training itself being carried out according to all applicable safety regulations.

### 1.2.3 Symbols used in the text and their meanings

**Warnings**  
are marked by a warning triangle

**Notes,**  
i.e. important information for the user,  
are marked by an information sign

(5)  
(Fig. 5) **Figures in brackets refer to item numbers in drawings / illustrations or to the drawings / illustrations themselves.**

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### 1.2.4 Liability

This document is strictly for the use of qualified service engineers with the requisite technical skills.

Only persons who have successfully completed the appropriate service training provided by Leica Microsystems GmbH, Nussloch and are in the employ of a company in the Leica Group or of an agency, distributor, or service workshop duly authorized by Leica Microsystems GmbH, Nussloch have the status of qualified service engineer.

Leica Microsystems GmbH, Nussloch accepts no liability whatever for direct or indirect damage that may occur due to the unauthorized or improper use or interpretation of this document by any person who is not a qualified service engineer in accordance with the above definition.

Service technicians have the following obligations:

- To understand and follow the safety information and instructions on the product and in the user manual.
- To be familiar with local regulations relating to industrial and non-industrial accident prevention in the knowledge that these regulations are up to date.
- To inform Leica immediately in writing if the equipment becomes unsafe.

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### 1.2.5 Instructions

To ensure trouble-free operation of the instrument at all times, the following instructions and warnings should be observed:

- The protective devices on the instrument and its accessories must not be removed or modified.
- Only service engineers authorized by Leica may access, service and repair the internal components of the instrument.

### 1.2.6 Instrument type

All information in this service manual applies only to the instrument type indicated on the title page. A nameplate with the serial number is fixed on the back of the instrument.

### 1.2.7 Transport and installation

- Do not operate the instrument in rooms where risk of explosion exists!
- Do not expose the instrument to direct sunlight (windows)!
- Do not install the instrument above a heater!
- Install the instrument on an even laboratory bench which must be absolutely level!
- Two people are needed to lift / carry the instrument!
- Before connecting the instrument to mains, make sure the correct voltage setting, matching the nominal voltage at the site of installation, has been selected!
- Upon installing the drain hose make sure there is a gradient from the drain outlet to the waste pipe.
- To protect the user from hazardous solvent fumes, make sure to operate the instrument either with the activated carbon filter or the exhaust air hose!

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### 1.3.1 Information on instrument design and safe handling

This instrument has been built and tested in accordance with the following safety regulations on electrical measuring, control, regulating and laboratory devices:

- DIN EN 292
- DIN EN 61010-1
- EN 50082-1
- EN 55011
- IEC 1000-4

as well as according to the international quality standard

- DIN ISO 9001

In order to maintain this condition and to ensure safe operation, the operator must observe the instructions and warnings contained in this instruction manual.

### 1.3.2 Integrated safety devices

The instrument is equipped with the following safety devices:

- Handwheel lock
- Handwheel grip centering
- Emergency stop function  
(instruments with sectioning motor only)
- Knife holder equipped with knife guard

The safety devices installed by the manufacturer of the instrument only constitute the basis of accident prevention.

Mainly responsible for accident-free operation is above all the institution which owns the instrument and, in addition, the designated personnel who operate, service or repair the instrument.

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### 1.3.3 Locking the handwheel


Always cover the cutting edge with the knife guard and lock the handwheel:

- Prior to doing any work on knife and / or specimen.
- Prior to exchanging specimens.
- During work breaks.

The handwheel can be locked in 2 positions:

- with the grip in the uppermost position (left),
- with the grip in the lowest position (right).

#### Locking

- Rotate handwheel, until grip (1) is in upper or lower position.
- To lock, press pin (2) to the right into position (2b). The upper locking position for pin (2) is marked by a black dot (4).

#### Instruments with sectioning motor

The sectioning motor is now blocked.

#### All instruments

The message 'LOCKED' in the display of control panel 1 indicates that the handwheel has been locked:

|   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |   |
|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|---|
| C | T | - | 3 | 0 | ° | C |   |   | O | T | - | 3 | 5 | ° | C |
|   |   |   |   |   |   | L | O | C | K | E | D |   |   |   |   |

#### Unlocking

- To unlock, push locking pin (2) to the left into position (2a).
- Display indication 'LOCKED' disappears.

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
## Instruments with sectioning motor

The sectioning motor can now be activated again.

### 1.3.4 Centering the handwheel grip


**During motorized sectioning, for safety reasons always center the handwheel grip!**


- To center grip (1), pull outwards and pivot into center of handwheel.
- When released, the grip locks into position.

### 1.3.5 Emergency stop function

The emergency stop is activated via the red emergency stop button in control panel 2 or via the footswitch

Control panel 2

Footswitch

#### Activating the emergency stop function

- Press emergency stop button or step on footswitch forcefully.  
STOP (red) lights up.  
As soon as the emergency stop function is activated, the sectioning motor stops.

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

**Deactivating the emergency stop**

- To deactivate, rotate emergency stop button in direction of arrow.

**If the emergency stop function has been activated by the footswitch, unlocking is not necessary (function is unlocked as soon as the footswitch is released).**

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**1.3.6 Knife guard**

All knife holders are equipped with a knife guard (see separate instruction manuals on knife holders).

**Always cover the cutting edge with the knife guard:**

- Prior to doing any work on knife and/ or specimen.
- Prior to exchanging specimens.
- During work breaks.

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**1.3.7 Transport**

- To avoid severe damage to the instrument by running it while the compressor oil is displaced from its regular position:
- Do not tilt the instrument, only transport in an upright position.
- After transport, wait at least 4 hours before turning the instrument on.  
The compressor oil may have been displaced during transport and must settle to its original position before switching the instrument on. Otherwise, the instrument may be severely damaged.

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**1.3.8 Site requirements**

- Do not operate the instrument in rooms with explosion hazard.
- To ensure proper instrument function maintain a minimum distance of 10 cm between walls and/ or furniture and all sides of the instrument!

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### 1.3.9 Electrical connections

- Do not use extension cords for connecting the instrument to mains.  
Fire hazard!  
Instrument malfunctions caused by voltage drop.
- During the start-up phase of the compressor, the nominal voltage must not drop below the values specified in chapter 1.3.31 'Technical data'!  
The compressor needs a start-up current at between 25 and 35 A .
- Ensure uniform current supply according to specifications. Electrical power supply deviating from specifications damages the instrument.
- Therefore, arrange for the electrical installations on site to be checked by a trained professional and make sure any necessary upgrades are installed!
- Have the circuit protected by a fuse of its own!
- Do not connect any other consumers to the same circuit.
- Prior to connecting the instrument to mains, make sure the electrical power supply in your laboratory corresponds to the values indicated on the instrument nameplate (located at the rear of the instrument)!

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### 1.3.10 Handling microtome knives/blades.

**Microtome knives and disposable blades have extremely sharp cutting edges and can cause serious injuries.**

- Handle knives / blades with utmost care.
- Never leave any knives / blades in unprotected places!
- Never place a knife, no matter where, with the cutting edge facing upwards.
- Never try to catch a falling knife.
- Always insert the specimen before inserting the knife.

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### 1.3.12 Motorized sectioning

- Do not interrupt sectioning / trimming by setting the sliding potentiometer to zero speed!  
This does not really switch the sectioning function off - it only operates at '0' speed. If the sliding potentiometer is accidentally moved, the instrument will resume sectioning immediately (risk of injury)!
- During motorized sectioning, always center the handwheel grip!

### 1.3.13 Defrosting/Handling frozen tissue

- Never leave specimens unattended in the cryochamber over an extended period of time! In case of power failure or instrument failure, or during the automatic defrost cycle, tissue material can be destroyed.
- During the defrost cycle the cryochamber is partially warmed. - Therefore:
- Remove sensitive specimens from the chamber prior to defrosting.
- If automatic defrosting is programmed to take place during the night, remember to remove all specimens from the cryo-chamber prior to leaving work.

### 1.3.14 Frozen parts of the instrument and frozen accessories

Prolonged contact of bare skin to frozen surfaces of the instrument or to frozen accessories (specimen discs, knife holder, shelves etc.) can cause frostbite.

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### 1.3.15 Infectious/radioactive material

**Use caution when working with potentially infectious specimens - Risk of infection!**

When working with potentially infectious / radioactive specimens:

- Wear protective clothes (gloves, protective boots, mask, lab coat), in compliance with radiation safety regulations and/or in-house regulations on handling infectious / radioactive material.

When working with radioactive specimens:

- Comply with applicable radiation safety regulations!
- Dispose of radioactive specimen waste according to applicable regulations.

### 1.3.16 Disinfection and cleaning

- Prior to disinfection, switch the instrument off and unplug it from mains.
- For removal of the microtome from the cryochamber, see chapter 1.3.17 Removing microtome.
- For disinfection, wear protective gear: (gloves, mask, lab coat etc.)!
- For disinfection, only use alcohol-based disinfectants!
- Do not use solvents (xylene, acetone etc.) for cleaning or disinfection!
- Do not spray disinfectants into the evaporator! Risk of icing!
- Explosion hazard when working with alcohol: Make sure the premises are appropriately ventilated!
- When using disinfectants and detergents, comply with all safety instructions supplied by the manufacturer of the product!
- Dispose of waste liquids from disinfection/cleaning as well as of sectioning waste according to applicable regulations on disposal of special category waste!
- Disinfected accessories must be thoroughly dry when reinserting them into the chamber! Risk of icing!
- Make sure the chamber is completely dry before switching the instrument back on: Explosion hazard through alcohol vapors! Risk of icing!

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### 1.3.17 Removing/reinstalling the microtome

#### Before removing the microtome

- Switch instrument off.
- Unplug from mains.
- Place handwheel grip in lowest position and lock.  
When removing the microtome, the specimen head must always be locked in the lowest position.  
Otherwise the upper part of the slot cover might be bent and consequently damaged!

#### When removing the microtome

- Wear gloves when removing the microtome while it is still frozen. Risk of frost bite!
- On instruments with specimen cooling: do not distort the refrigerating tube!  
If distorted it might break, causing extremely cold refrigerant to escape. Risk of frost bite!

#### Before reinstalling the microtome

- Microtome must be completely dry.  
Humidity in the interior of the microtome freezes and causes microtome malfunctions and/or damage to the microtome.
- All accessories/tools removed from the cryochamber must be thoroughly dry before putting them back into the chamber!  
Risk of icing!

### 1.3.18 Display message 'Dry microtome'

If the error message 'Dry Microtome' is displayed in control panel 1, the following has happened:

- Cryochamber refrigeration has been interrupted for an extended period of time (e.g. power failure), causing the chamber temperature to rise into the positive digits.
- If this message appears, do not switch on the instrument but remove the microtome from the chamber, disinfect, if necessary, and dry thoroughly before reinstalling it into the chamber.

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### 1.3.19 Maintenance

Only technical service engineers authorized by Leica may access the internal components of the instrument for service and repair. The fluorescent light lamp (chamber illumination), unless broken or splintered, can be replaced by the user:

- Switch off mains switch!
- Unplug the instrument from mains!
- If the lamp is broken or splintered: Have lamp replaced by Leica Technical Service! Risk of injury!
- Use only those replacement lamps that correspond to technical specification (see chapter 1.3.31 'Technical Data').

### 1.3.20 Site requirements

**Make sure to read and follow all safety instructions in chapter 1.1.8 'Site requirements'!**

### 1.3.21 General site requirements

- No direct sunlight.
- Electrical power supply within distance (length of power cord = approx. 4 meters - do not use extension cords!).
- No draft (caused by air conditioning etc.).
- Even floor surface.
- Practically vibration-free floor.
- Handwheel easily accessible.
- Room temperature constantly below +22 °C.
- Relative humidity of air maximum 60 %.

**High ambient temperature and/or high air humidity negatively affect instrument cooling performance!**

### 1.3.22 Electrical connections

**Make sure to read and follow all safety instructions in chapter 1.3.9 'Electrical connections'!**

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### 1.3.23 Unpacking and installation

Unpacking instructions are always located in a transparent protective envelope on the outside of the instrument shipping crate.

**Make sure to read and follow all safety instructions provided in chapter 1.3.7 'Transport' and on the unpacking instructions!**

### 1.3.24 Repacking

We recommend to keep the original shipping crate and the unpacking instructions for the Leica CM3050 S. For repacking, proceed as per unpacking instructions, in reverse order.

### 1.3.25 Standard delivery

Together with the Leica CM3050 S instrument you receive a box containing the following standard accessories:

- 1 Handwheel (including screw, spring washer and cover disc)
- 1 Heat extractor, stationary
- 1 Low-temperature stabilizer for heat extractor
- 3 Specimen discs, 30 mm ø
- 1 Storage shelf, right
- 1 Storage shelf, left
- 1 Section waste tray
- 1 Rubber mat
- 1 Cover for quick-freeze shelf
- 2 Brushes (brush, fine and 'Leica' brush)
- 1 Shelf for brushes
- 1 Tool set
  - 1 Allen key, size 4
  - 1 Allen key, size 5
  - 1 Allen key, size 6
  - Single-ended open-jaw wrench
- 1 Bottle of embedding medium for cryosectioning, 125 ml
- 1 Bottle of cryostat oil No. 407, 100 ml
- 1 Instruction manual  
(German, English, French, Spanish)

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**In addition to the above**

Instruments with specimen cooling:

- 1 90° Prism  
for low temperature quick freezing

Instruments with sectioning motor:

- 1 Footswitch with protective guard

Configured instruments:

- 1 Knife holder base
- 1 Knife holder with accessories

**Further accessories**

Further accessories which you ordered will be included in the box containing the standard delivery items.

**Knife holders are delivered with anti-roll guide, knife guard, and a separate instruction manual.**

**In case of non-configured instruments, the knife holder is not a part of standard delivery but must be ordered separately.**

Check all delivered parts against the packing list and against your order to verify whether the delivery is complete! If there is any difference, contact your local Leica office immediately.

**Available models**

- Basic instrument  
with sectioning motor  
without specimen cooling
- Basic instrument  
without sectioning motor  
with specimen cooling
- Basic instrument  
with sectioning motor  
with specimen cooling

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
### 1.3.26 Installing the handwheel

- Insert pin (1) of the handwheel shaft into hole (2).
- Place spring washer (3) onto screw (4) as shown.
- Tighten screw (4) with an Allen key.
- Cover with selfadhesive disc (5).

**For purposes of transport (e.g. narrow doors), the handwheel can be removed.**

To remove the handwheel, proceed as described above but in reverse order.

### 1.3.27 Inserting the accessories


- Place the rubber mat on top of the housing.
- Insert the storage shelves into the cryochamber.
- Install the stationary heat extractor into the quick-freeze shelf.
- Insert the low temperature stabilizer into the quick freeze shelf (it must be located in the pivoting range of the heat extractor).
- Insert section waste tray and brush shelf.
- Install knife holder base onto microtome base plate and clamp.
- Install knife holder and clamp (see knife holder instruction manual for details).
- Place knife case with knife into chamber to precool.
- Place all tools needed for section preparation into the chamber.
- Close the sliding window.

**For a complete overview of all individual parts, see Chapter 1.3.30 'Overview'.**

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
### 1.1.28 The footswitch

#### Function

The footswitch performs the same task as the RUN/STOP and RUN/ENABLE keys (activating/deactivating motorized sectioning / trimming). In addition, the footswitch can be used to activate the emergency stop function.

#### Models with footswitch


All instruments with sectioning motor.


**In all instrument models that are delivered with footswitch, the footswitch must be installed! - Otherwise the instruments are not functional.**

#### Connecting the footswitch

- Insert footswitch into port (1) and secure.

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
**1.3.29 Prior to switching on the instrument**

**After transport, observe a waiting period of at least 4 hours before turning the instrument on! - See also safety instructions in chapter 1.3.7 'Transport'.**

**Have you observed all safety instructions in chapters 1.3.8 'Site requirements' and 1.3.9 'Electrical connections'? If not, Please read chapters 1.3.8 and 1.3.9!**

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


## 1.3.30 Overview


- 1 Cryostat CM3050 S
- 2 Control panel 1
- 3 Control panel 2
- 4 Storage shelf, left
- 5 Storage shelf, right
- 6 Rubber mat
- 7 Mains switch /  
Automatic cutout for instrument,  
Automatic cutout for sectioning  
motor,  
Footswitch port
- 8 Quick-freeze shelf
- 9 Stationary  
heat extractor
- 10 Mobile heat extractor
- 11 Specimen disc
- 12 Thermoblock (optional)
- 13 Section waste tray
- 14 Brush shelf
- 15 Specimen head w/o specimen  
cooling
- 16 Specimen head with  
specimen cooling (Option)
- 17 90° Prism  
(instruments with specimen cooling only)
- 18 Knife holder base
- 19 Knife holder CE
- 20 Knife holder CN
- 21 Knife holder CS
- 22 Footswitch with  
protective guard

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## 1.3.31 Technical Data

**Operating temperature range (ambient temperature): + 18°C to + 40°C. All specifications related to temperature are valid only up to an ambient temperature of + 22°C and relative air humidity of less than 60%!**

| Type                              | CM3050S-10          | CM3050S-1           | CM3050S-8           | CM3050S-9           |
|-----------------------------------|---------------------|---------------------|---------------------|---------------------|
| Marks of conformity               | -                   | -                   | -                   | -                   |
| Nominal voltage                   | 100 V AC $\pm 10\%$ | 230 V AC $\pm 10\%$ | 240 V AC $\pm 10\%$ | 100 V AC $\pm 10\%$ |
| Nominal frequency                 | 50 Hz               | 50 Hz               | 50 Hz               | 60 Hz               |
| Power draw                        | 1800 VA             | 1800 VA             | 1800 VA             | 1800 VA             |
| Max. start-up current for 5 sec   | 35 A eff.           | 25 A eff.           | 25 A eff.           | 30 A eff.           |
| Protective fuse                   | I                   | I                   | I                   | I                   |
| Automatic cutout                  | T15A M3             | T10A T1             | T10A T1             | T15A M3             |
| Pollution degree ①                | 2                   | 2                   | 2                   | 2                   |
| Overvoltage installation category | II                  | II                  | II                  | II                  |
| Heat emission (max.)              | 1800 J/s            | 1800 J/s            | 1800 J/s            | 1800 J/s            |

| Refrigeration system             | 50 Hz                                                                                                                                | 60 Hz                                                                                                                                 |
|----------------------------------|--------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| <b>Cryochamber</b>               |                                                                                                                                      |                                                                                                                                       |
| Temperature selection range      | 0° C to -40°C (+ 2 K / - 0 K)                                                                                                        | 0° C to -40°C (+2 K / - 0 K)                                                                                                          |
| Defrosting                       | Automatic hot-gas defrost cycle, programmable in 15-min increments; 1 automatic defrost cycle per 24 hours. Duration: 6 - 12 minutes | Automatic hot-gas defrost cycle, programmable in 15-min. increments; 1 automatic defrost cycle per 24 hours. Duration: 6 - 12 minutes |
| Refrigeration capacity ②         | Manual defrost cycle<br>690 W                                                                                                        | Manual defrost cycle<br>690 W                                                                                                         |
| Safety factor                    | 3                                                                                                                                    | 3                                                                                                                                     |
| Refrigerant                      | 280 g $\pm 5$ g refrigerant R 404A *                                                                                                 | 280 g $\pm 5$ g refrigerant R 404A *                                                                                                  |
| Compressor oil                   | 0.6 l EMKARATE RL-22S, ICI *                                                                                                         | 0.6 l EMKARATE RL-22S, ICI *                                                                                                          |
| <b>Quick-freeze shelf</b>        |                                                                                                                                      |                                                                                                                                       |
| Maximum temperature              | - 43°C (+ 0 K / - 2 K)                                                                                                               | - 43°C (+ 0 K / - 2 K)                                                                                                                |
| Number of quick-freeze stations  | 10                                                                                                                                   | 10                                                                                                                                    |
| <b>Specimen cooling (Option)</b> |                                                                                                                                      |                                                                                                                                       |
| Temperature range                | -10° C to -50°C $\pm 3$ K                                                                                                            | -10° C to -50°C $\pm 3$ K                                                                                                             |
| Defrosting                       | Manual defrost cycle (electric)                                                                                                      | Manual defrost cycle (electric)                                                                                                       |
| Refrigeration capacity ②         | 320 W                                                                                                                                | 320 W                                                                                                                                 |
| Safety factor                    | 3                                                                                                                                    | 3                                                                                                                                     |
| Refrigerant                      | 185 g $\pm 5$ g refrigerant R 404A *                                                                                                 | 185 g $\pm 5$ g refrigerant R 404A *                                                                                                  |
| Compressor oil                   | 0.4 l alpha 22, Kyodo *<br>or RENISO E22, Fuchs *                                                                                    | 0.4 l alpha 22, Kyodo *<br>or RENISO E22, Fuchs *                                                                                     |

**\*) Refrigerant and compressor oil to be replaced only by trained and authorized service personnel!**

|                                                                                                                                                                                                                                                                                                                                                                                                   |                                                                                                                                                                                                     |                                                                                                                                                                                                     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <div><div>CM3050S-3</div><div>c-UL</div><div>120 V AC ±10%</div><div>60 Hz</div><div>1800 VA</div><div>35 A eff.</div><div>I</div><div>T15A T1</div><div>2</div><div>II</div><div>1800 J/s</div></div>                                                                                                                                                                                            | <div><div>CM3050S-6</div><div>-</div><div>208 V AC ±10%</div><div>60 Hz</div><div>1800 VA</div><div>30 A eff.</div><div>I</div><div>T12A T1</div><div>2</div><div>II</div><div>1800 J/s</div></div> | <div><div>CM3050S-7</div><div>-</div><div>230 V AC ±10%</div><div>60 Hz</div><div>1800 VA</div><div>25 A eff.</div><div>I</div><div>T10A T1</div><div>2</div><div>II</div><div>1800 J/s</div></div> | <div><div>Microtome</div><div>Rotary microtome</div><div>Section thickness setting</div><div>Overall specimen feed</div><div>Vertical stroke</div><div>Specimen retraction</div><div>Maximum specimen size</div><div>Specimen orientation</div><div>Trimming (6 thickness settings)</div><div>Motorized coarse feed</div><div>- slow</div><div>- fast</div><div>0.5 - 300 µm</div><div>25 mm</div><div>59 mm</div><div>50 µm</div><div>40 x 55 mm</div><div>8° (x-, y-, z-axis)</div><div>5 - 150 µm</div><div>500 µm/s</div><div>1,000 µm/s</div></div> |
| <div><div>Fluorescent light lamp (cryochamber illumination):</div><div>50 Hz version:</div><div>Osram Dulux S 11 W/21</div><div>Light color 'LUMILUX hellweiß'</div><div>(brilliant white)</div><div>60 Hz version:</div><div>Osram Dulux S 13 W/21</div><div>Light color 'LUMILUX hellweiß'</div><div>(brilliant white)</div></div>                                                              |                                                                                                                                                                                                     |                                                                                                                                                                                                     | <div><div>Sectioning motor (Option):</div><div>Sectioning speed</div><div>- slow</div><div>- fast</div><div>min.</div><div>max.</div><div>v<sub>max</sub></div><div>min.</div><div>max.</div><div>v<sub>max</sub></div><div>0.1 mm/s</div><div>100 mm/s</div><div>210 mm/s</div><div>0.1 mm/s</div><div>170 mm/s</div><div>210 mm/s</div></div>                                                                                                                                                                                                          |
| <div><div>Cryocabinet</div><div>Dimensions</div><div>Width (w/o handwheel)</div><div>Width (including handwheel)</div><div>Depth</div><div>Height (arm rest level)</div><div>Overall height</div><div>800 mm</div><div>882 mm</div><div>766 mm</div><div>840 mm</div><div>1040 mm</div></div> <div><div>Weight</div><div>(incl. microtome, w/o specimen cooling)</div><div>ca. 180 kg</div></div> |                                                                                                                                                                                                     |                                                                                                                                                                                                     | <div><div>①</div><div>according to IEC-1010, UL 3101</div><div>②</div><div>according to CECOMAF:</div><div>Liquid temperature 45°C</div><div>Evaporation temperature -25°C</div></div>                                                                                                                                                                                                                                                                                                                                                                   |

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

**Disassembly - control panel, small (2)**

**Before performing any work on electrical components, unplug the instrument!**

1. Disconnect all cable connections of control panel, small (2), see wiring diagram.
2. Remove nuts (12) together with tooth lock washers (11) and washers (13).
3. Remove the small keyboard (2) upwards.

**Disassembly - control panel, large (3)**

**Before performing any work on electrical components, unplug the instrument!**

1. Disconnect all cable connections of control panel, large (3), see wiring diagram.
2. Remove nuts (12) together with tooth lock washers (11) and washers (13).
3. Remove control panel, large (3) upwards.

**Disassembly  
- cover plate, mains switch unit (4)**

**Before performing any work on electrical components, unplug the instrument!**

1. Unscrew screws (21) and remove cover (20).
2. Disconnect all cable connections of the mains switch unit cover plate (4), see wiring diagram.
3. Remove nuts (12) together with tooth lock washers (11) and washers (25).
4. Remove the mains switch unit cover plate (4) forward.
5. Remove screws (22) and safety cutout (10).
6. Remove nut (9.2) and washer (9.1). Next, remove safety cutout (9).
7. Remove distance bolt (24) with washers (13) and take off filter board (23).

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
Abb. 090999k2

**Disassembly - drain pan (19)**

1. Remove screw (22) together with retaining washer (21) and retaining ribbon (20).
2. Remove the drain pan (19).

**Disassembly - guard plate (1)**

**Before performing any work on electrical components, unplug the instrument!**

1. Disconnect all cable connections located at the guard plate (1), see wiring diagram.
2. Remove the guard plate (1) upwards.

**Disassembly - leaf spring (12)**

1. Remove screw (14) and tooth lock washer (13).
2. Take off leaf spring (12).

**Disassembly - reflector (7)**

**Before performing any work on electrical components, unplug the instrument!**

1. Disconnect all cable connections at socket (9), see wiring diagram.
2. Remove fluorescent lamp (16) to the left.
3. Unscrew screws (3) and remove the reflector (7) upwards.
4. Remove screws (4) and washers (5) as well as insulating angle (10).
5. Remove socket (9), holder (17) and cable clip (8) downwards.

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
Abb. 090999k3

**Disassembly - microtome with base plate (7)**

1. Slightly lift the heated window (13) and pull it out towards the front.
2. Pull off sensor (5) backwards.  
**Do not pull at the cable!**
3. Take the coarse feed cable (6) at the plug and pull it to the left. **Do not pull at the cable!**
4. On instruments with specimen cooling:
  - a) Remove screw (1).
  - b) Slightly pull specimen head (4) forward.
  - c) Use the specimen head fixture (11) to hang up the specimen head (4) on the drain pan (12).
5. Remove screws (2) and saucer springs (3).
6. Turn the handwheel (14) to move the specimen cylinder (15) to its lower end position.
7. Move the microtome with base plate (7) to the left as far as it will go.
8. Remove microtome / base plate (7) forward.

**Reinstallation - microtome / base plate (7)**

1. Reinstallation is done in reverse order.
2. To reconnect the components of the claw coupling (8) the specimen cylinder (15) and the handle of handwheel (14) must be each at their lower end position.

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
Abb. 090999k5

**Disassembly - housing (1)**

1. Remove screw (3) and pull housing (1) forward far enough so that slot cover (2) can be removed upwards.
2. Take off housing (1).

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**Disassembly - sensor clip (6)**

**Before performing any work on electrical components, unplug the instrument!**

1. Pull off sensor (9) backwards.  
**Do not pull at the cable!**
2. Remove screw (4) and washer (5).
3. Remove sensor clip (6) and spring (11).

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**Disassembly - sensor (9)**

**Before performing any work on electrical components, unplug the instrument!**

1. Disconnect all cable connections at sensor (9), see wiring diagram.
2. Pull off sensor (9) backwards.  
**Do not pull at the cable!**

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**Disassembly - coarse feed cable (10)**

**Before performing any work on electrical components, unplug the instrument!**

1. Disconnect all cable connections at coarse feed cable (10), see wiring diagram.
2. Take coarse feed cable (10) by the jack plug and pull it off backwards. **Do not pull at the cable!**

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
Abb. 090999k6

**Disassembly - eccentric bolt (13)**

1. Remove set screw (14).
2. Rotate clamping lever (16) either backward or forward while at the same time pulling it out sideways.
3. Remove clamping sleeve (15) and then clamping lever (16).

**Disassembly - pulling bolt (11)**

1. Remove screw (18) then remove clamping device (23) and saucer springs (19).
2. Disassemble eccentric bolt (13), see above.
3. Pull out pulling bolt (11) upwards.
4. Pull off pressure spring (8).

**Disassembling angle steel, short (5)**

1. Detach screws (27) and remove washers (28).
2. Detach screws (3), remove angle steel, short (5) and washers (6).

**Disassembly - angle steel (25)**

1. Loosen screws (27) and remove washers (28).
2. Remove angle steel (25) forward.
3. Remove screws (27) and detach washers (28) together with angle steel (26).
4. To adjust the microtome centrically to the motor drive, move angle (25) up and down and angle (26) back and forth. Tighten screw (22) just slightly - it serves as a limit point.

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
**Disassembly - cover plate, right (12) and cover plate, left (13)**

1. Unhook tensions springs (14).
2. Detach screws (6) and remove spring suspensions (16).
3. Remove spring hooks (15).
4. Detach screws (3) and remove cover plate, right (12) as well as cover plate, left (13).

**Disassembly - slider (2) and cross roller guides (10)**

1. Remove the housing (see chapter 2.4.3, item no.1).
2. Remove specimen head (see chapter 2. 4.1, item no. 4), slot cover (see chapter 2.4.3, item no. 2), micrometer mechanism, assy. (see chapter 2.4.6) and cover plates (13) and (12).
3. Place frame (1) with base plate (see chapter 2.4, item no. 4) on its rear side. Support rear of frame (1) with one hand while tilting it.
4. Detach setscrews (5) and hexagon nut (8).
5. Detach screws (4) which hold cross roller guide rail (I) and remove slider (2) forward.
6. Next, remove screws (4), which hold cross roller guide rails (II), (III) and (IV), from slider (2) and frame (1).

**Adjusting slider (2) and cross roller guides (10)**

1. Press cross roller guide rails (II) and (III) with your fingers into the square angles at the front side of the slider (2) and tighten screws (4). Torque = 150 Ncm.
2. Press cross roller guide rail (IV) into the square angle at the front of frame (1) and fasten with screws (4). - Torque = 150 Ncm.
3. Place roller cages (V) in position.

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

Abb. 250899k1

**Disassembly - micrometer mechanism**

1. Remove the microtome (see chapter 2.4.1).
2. Remove the housing (see chapter 2.4.3).
3. Remove six screws (26) from the front of slider (21), pull out the entire micrometer mechanism backwards.
4. Unhook spring (19).
5. Loosen clamping screw (12) from the flexible membrane coupling which clamps the stepper motor shaft.
6. Remove screws (31) and washers (32), which hold the microswitches, as well as screws (38) and washers (39), which hold cable clamps (42 and 40).
7. Remove screws (35), washers (30) and spacer (11). Then detach the feed mechanism assembly (2).
8. Loosen the clamping screw of flexible membrane coupling (12), which clamp spindle shank (1.1).
9. Remove flexible membrane coupling (12) from spindle shank (1.1).
10. Loosen screw (37) and remove clamping ring (13); saucer spring (14) and deep groove ball bearing (15).
11. Remove screws (29) and washer (30); detach bearing washer (5) from spindle shank (1.1).
12. Pull off deep groove ball bearing (15).
13. Unscrew (9) V-ring.
14. Pull spindle (1.1) - together with nut (1.2), needle bearing (16) and bushing (8) - out of cylinder (3).
15. Remove screw (33) and washer (34); then remove shifting block (6).
16. Remove screws (27) and washers (28) and pull ellipse (22) out of cylinder guide (4).
17. Pull cylinder (3) out of cylinder guide (4).
18. Remove screws (29) and washers (30); then remove bar (24).

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
Abb. 250899k1

**Adjustment - sliding block (41)**

1. Place bar (24) in its correct position; add screws (29) and washers (30) - tighten just slightly.
2. Place sliding blocks (41) in their guide rail.

**Attention! Make sure the sliding blocks are positioned correctly. The wide facet must face downward and be located on the right; the mark on sliding block (41) must be on the top right.**

- 3.) Move bar (24) up and down to adjust the coefficient of friction.

**At an ambient temperature of +23 °C the coefficient of friction should be about 2600 g. The coefficient of friction is temperature-dependent.**

- 4.) Tighten screws (29). - Torque = 350 Ncm.

**Assembly - micrometer mechanism**

1. Clean cylinder (3) and cylinder guide (4).
2. Insert spindle (1.1) and nut (1.2) into cylinder (3) and push them inward as far as they will go.

**Attention! Checking the clearance of the spindle. This check can be done only with a completely assembled microtome. The check should be performed at a cryochamber temperature of approx. -25 °C.**

- Fasten the dial gauge of 0.001mm scale division, part. no. 0178 11205, with its holder, part no. 0800 35114, to slider (21).
- Establish contact between the dial gauge and cylinder (3).
- Press a force measuring device from the front against cylinder (3) (with a force of 5000 g) and hold that value.
- Log the measured value into the dial gauge.
- Now use the force measuring device to pull cylinder (3) forward at a value of 5000 g and hold.
- Log that value into the dial gauge as well.
- The difference between the two logged values must be in the range from 5µm to 35µm
- This check has to be repeated at several different spindle / nut positions.

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
Abb. 250899k1

19. Assemble spacers (11) and fasten screws (35) and washers (30) just slightly.
20. Center feed mechanism assembly (2) (shift until centered) and tighten screws (35) crosswise. Torque = 200 Ncm.
21. Tighten the screws of the flexible membrane coupling.
22. Install both microswitches (2.2) parallel to the cylinder guide (4) and fasten with screws (31) and washers (32). - Torque =  $11 \pm 1$  Ncm.
23. Put cable clamps (40 and 42) in correct position and fasten with screws (38) and washers (39).
24. Hook in tension spring (19).
25. Measure torque IV.

**Attention! Torque IV is temperature-dependent and may be max. (see table below).**

+ 23°C -> 10Ncm

- 25°C -> 12Ncm

- 35°C -> 15Ncm

It is absolutely necessary to verify these values, as values exceeding those listed above may lead to problems in the feed mechanism.

26. Insert the micrometer mechanism assembly from the back into the opening in slider (21) and fasten it with screws (26).

**Attention! Make sure the micrometer mechanism is exactly centered (adjust if necessary)!**

27. Install the microtome in the chamber and cool down.
28. Do final check.

**Attention! The final check must be performed at a cryochamber temperature of approx. -25°C and must include the following steps:**

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
Abb. 250899k1

- Check the feed mechanism. For that purpose, mount dial gauge, part no. 0178 11205 and dial gauge holder, part no. 0800 35114. Subsequently (via the electronics), set the feed thickness values listed below. Perform 10 handwheel rotations at each setting. After each sequence of 10 rotations, the following feed values have to be measured:

|      |    |    |   |       |
|------|----|----|---|-------|
| 1µm  | -> | 14 | - | 16µm  |
| 2µm  | -> | 23 | - | 27µm  |
| 5µm  | -> | 45 | - | 55µm  |
| 10µm | -> | 90 | - | 110µm |

When doing these measurements, the specimen cooling must be "OFF".

- Retraction of the micrometer mechanism has to be  $50\mu\text{m} \pm 20\mu\text{m}$ .
- Coarse feed travel time 'Fast forward/backward' from 'Home' to 'Stop' has to be 20 - 24sec.

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
Abb. 141200k1

**Disassembly - shaft, assy. (6)**

1. Detach screws (14) and remove clutch (15) together with clutch components (15.1) and (15.2).
2. Remove dowel pin (13) and then ring (18) together with washer (7).
3. Remove screws (11) with washers (10).
4. Pull off bearing flange (9). Slide bearing (12) is friction-set into bearing flange (9).
5. Pull off deep groove ball bearing (8) and washer (7).
6. Pull out shaft, assy. (6) sideways, then remove pin (17) and sliding block (5).
7. See chapter 2.4.6 on how to adjust sliding block (5).
8. To reassemble, proceed in reverse order.

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
1. Disconnect the cable connections at microswitch (26) and precision potentiometer (48), see chapter 'Electronics'.
2. Remove screws (2) and retaining washers (3) from handwheel assy. mounting plate (1).
3. Remove hexagon nuts (4) and washers (3) from mounting plate (1).
4. Unhook toothed belt (52) and detach mounting plate (1).

**Disassembly - precision potentiometer (48)**

1. Loosen setscrew (47) and remove spur wheel (46).
2. Remove screws (33) and washers (45).
3. Remove precision potentiometer (48) from potentiometer elbow (43) .

**Disassembly - bearing flange (32)**

1. Remove locking washer (11) from primary shaft (38).
2. Detach setscrew (44).
3. Remove toothed wheel (42) and tooth lock washer Z5.1 (40).
4. Remove locking washer (37) and pull off primary shaft (38) forward.
5. Remove bearing flange (32) from bearing plate (31).
6. Remove both deep groove ball bearings (35) together with levelling washer (36).

**Keep in mind when reassembling:**

1. **Bearing flange (32) has to be glued into bearing plate (31) with Loctite 638.**
2. **In case that toothed wheel (42) had been removed from tooth lock washer (40), screws (41) have to be glued into tooth lock washer (40) with Loctite 638.**

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
**Disassembly - roller (54)**

1. Remove screw (56) and retaining washer (3).
2. Remove tenon block (58) and axle (57).
3. Remove deep groove ball bearing (55) from roller (54). - **Attention - secured with Loctite 638!**

**Disassembly - switch holder (21)**

1. Remove screws (24) and locking washers (23).
2. Remove switch holder (21) with components.
3. Unhook tension spring (28).
4. Detach screw (29) and remove switch lever (22).
5. Remove screws (27) and microswitch (26).
6. Remove grooved pin (25).

**Keep in mind when reassembling:**

1. Screw (29) has to be glued into switch holder (21) with Loctite 638.
2. Switch lever (22) has to be smooth-running but free from play.

**Disassembly - axle (6)**

1. Remove hexagon nut (19) and washer (20).
2. Detach axle (6) and attached components and remove forward.
3. Remove locking washers (14 and 11).
4. Remove washer (10), levelling washer (9) and levelling washer (15).
5. Pull off wheel hub (12) and crown gear (13) towards the front.

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
Abb. 960716r2

6. Remove both deep groove ball bearings (7) and then bush (8).

**Keep in mind when reassembling:**

1. In case crown gear (13) had been removed from wheel hub (12), crown gear (13) has to be glued onto wheel hub (12) with Loctite 638.
2. In case parallel pin (17) had been pulled out, it has to be glued into wheel hub (12) with Loctite 638.

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1. Disconnect the cable connection at the motor (5), see chapter 'Electronics'.
2. Disconnect cable connection at the bearing block (28), see chapter 'Electronics'.
3. Detach screw (see chapter 2.5, item no. 56).
4. Slightly push roller (see chapter 2.5, item no. 54) downwards and unhook toothed belt (44).
5. Remove nuts (2) and retaining washers (3).
6. Remove mounting plate (1) with all attached components.

**Disassembly - motor (5)**

1. Remove screws (7), locking washers (8) and distance sleeves (6).
3. Remove motor (5) towards the rear.
4. Take off toothed belt (12).
5. Remove setscrew (10) and then tooth lock washer Z1 (11).
6. Detach screws (64).
7. Pull off connector (65) and distance sleeve (66).

**Disassembly - bearing block (28)**

1. Take off locking washer (45).
2. Loosen setscrew (10) take off tooth lock washer Z5 (43).
3. Remove nuts (31) with retaining washers (3).
4. Unhook toothed belt (26).
5. Remove the entire bearing block with all attached components (28).
6. Loosen setscrew (52) and remove electromagnetic clutch (48).

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#### Remember when reassembling:

1. The intermediary part (48a) of the clutch may have a maximum axial play of 0.2 mm. If this is not the case, compensate with adjusting washer (50).
2. In case the threaded rods (30) had been loosened, they have to be glued into the bearing block (28) with Loctite 638.
3. Glue tooth lock washer Z4 (38) onto bearing bush (37) with Loctite 638.

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### Disassembly - roller (54)

1. Remove screw (57) and retaining washer (3).
2. Take out roller (54) with deep groove ball bearing (56) and axle (59).  
**Attention, deep groove ball bearing (56) is glued into roller (54) with Loctite 638.**
3. Pull ball bearing (56) / axle (59) out of roller (54).

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
Abb. 960717r2

**Disassembly - bearing flange (15)**

1. Remove locking washer (20) and then pull out to the rear intermediary axle (21) with tooth lock washer Z2 and Z3 (22 and 24).
2. Take off deep groove ball bearing (18) and remove levelling washer (19).
3. Remove screws (16) with washers (17) and detach bearing plate (14).
4. Pull off deep groove ball bearing (18).

**Keep in mind when reassembling:**

1. Glue bearing flange (15) into bearing plate (14) with Loctite 638.

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**Disassembly tooth lock washers (22 and 24)**

1. Remove locking washer (25) and parallel pin (23).
2. Detach tooth lock washers (22 and 24) from intermediary axle (21).

**Keep in mind when reassembling:**


1. Glue tooth lock washer (22) onto intermediary axle (21) with Loctite 638.

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**Detaching the handwheel**

1. Remove cover disc (19).
2. Unscrew screw (16) from axle (see chapter 2.5, item no. 6).
3. Remove handwheel (10).

**Disassembly - handwheel**

1. Remove hexagon nuts (7 and 8) with discs (9).
2. Detach plate (18) to the front, remove screws (27) and washer (20).
3. Take off sliding pin (12), washer (13), disc (17), screw (16) and conical spring washer (15).
4. Remove screws (14) and pull off handle, assy. (26).
5. Unscrew shaft (22), washer (23), pressure spring (24) and bearing (25) from handle (26).
6. Remove screws (11) and detach bearing (2).
7. Remove screw (6) with hexagon nut (4) and locking washer (3).
8. Take off leaf spring (5) and bolt (1).

**Remember when reassembling:**

1. Shaft (22) is glued into handle (26) with Loctite 274.

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

**Detaching the handwheel**

1. Remove cover disc (12).
2. Remove screw (11) and conical spring washer (10).
3. Take off handwheel (9).

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**Disassembly - handwheel**

1. Unscrew handle (15) from the back using an Allen key.
2. Remove screw (7) and locking washer (8).
3. Take off bearing (2).
4. Remove screw (6) with hexagon nut (4) and locking washer (3).
5. Remove leaf spring (5) and bolt (1).

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
Abb. 960718r2

### Disassembly - Specimen Orientation

1. Rotate screw (1) until specimen orientation, assy., can be removed from cover (19).
2. Pull out parallel pins (16) and remove spindles (12) and (14).
3. Remove nut (18) from universal joint (17).
4. Loosen setscrew (4).
5. With a size 3 Allen key unscrew clamping lever (23) from setscrew (21).
6. Remove washer (22).
7. Remove Ribe caps (2) and unscrew screws (3).
8. Remove ball socket (5) together with pressure springs (9), (8) and screws (1).
9. Remove ball (7) together with universal joint (17).  
**Universal joint (17) is glued into ball (7) with Loctite 638!**
10. Detach countersunk screws (20) and remove cover (19).

### Reassembling the specimen orientation

1. Place ball (7) together with universal joint (17) into ball socket (11).
2. Insert springs (9) and (8) and screws (1) into their correct position.
3. Insert setscrew (4) into ball socket (5) and tighten.
4. Place ball socket, front (5) onto springs (8) and (9) and fasten just slightly with screw (3).
5. Insert setscrew (21) from the back into the drilling in ball socket (11), add washer (22) and screw on clamping lever (23).  
**Rotate setscrew (4) and setscrew (21) as necessary to adjust ball socket (5) relative to ball socket (11). The two sockets have to be parallel to each other!**  
**Attention! - Nut (15) is secured with Loctite 638!**

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
Abb. 960718r2

6. Tighten screws (3).
7. Rotate setscrew (21) further to adjust the clamping of ball (7).  
Coefficient of friction: 2.5 kg  
Measuring point: outermost point of clamping lever (23).
8. Tighten setscrew (4).
9. Insert nut (18) into universal joint (17).
10. Insert spindles (14) and (12) with knurled knobs (13) into the drilling in ball socket (11) and secure with parallel pin (16).  
Make sure both spindles (12) and (14) are positioned correctly!  
Also, both parallel pins (16) must have some clearance within their guide bush!

**If the parallel pins (16) get jammed, spindles (12) and (14) can no longer be disassembled, neither can the clamping of ball socket, front (5) be adjusted any more!**

11. With screws (1) fasten specimen orientation, assy. to cover (19).

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
Abb. 051099k1

**Disassembly - heat extractor**

1. Remove screw (2) and then bolt (1). Screw (2) is secured with Loctite 274.
2. Remove screw (5); then remove button (6) in upward direction and pressure spring (8) together with axle (9) and stamp (10) in downward direction.

**Attention! Axle (9) is glued into stamp (10) with Loctite 638.**

3. Loosen screw (4) and pull bushing (7) out of clamp (3).
4. For reassembly proceed in reverse order.

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

(\*) only with sectioning motor - otherwise dummy plug RM 2145, part-no. 6851 01001, on periphery board

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



Abb. 044333557-3c

**General Information:**

The foot switch is identical to the one used in the RM21..series. It switches the cutting motor on/off and may also be used to activate the emergency stop.

The foot switch filter filters all disturbances, which could otherwise verge from the foot switch into the control electronics.

Replacement foot switches are always delivered with a gender changer!

When connecting the replacement foot switch to a CM3050 S, do not use the Gender Changer!

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
**General Information:**

The power supply unit filters the mains voltage and distributes it to the transformers. For some mains voltages a superpose transformer is needed for the compressor drive.

All compressors for the refrigerating system are chosen only according to mains frequency.

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50Hz      ⇒      230V Compressor  
60Hz      ⇒      115V Compressor

**Table: Assignment of transformers and compressors**

| Mains voltage ( V ) | Frequency ( Hz ) | Compressor Voltage ( V ) | Superpose transformer No. |
|---------------------|------------------|--------------------------|---------------------------|
| 230                 | 60               | 115                      | 0443 24875                |
| 208                 | 60               | 115                      | 0443 24875                |
| 200                 | 60               | 115                      | 0443 24875                |
| 127                 | 60               | 115                      | 0443 24875                |
| 120                 | 60               | 115                      | 0443 24875                |
| 115                 | 60               | 115                      | 0443 24875                |
|                     |                  |                          | 0443 24875                |
| 240                 | 50               | 230                      |                           |
| 230                 | 50               | 230                      | 0443 28640                |
| 200                 | 50               | 230                      | 0443 28640                |
|                     |                  |                          |                           |
| 100                 | 60               | 115                      | 0443 29091                |
| 100                 | 50               | 230                      | 0443 29103                |

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
Leica Instruments  
699000057  
92-01-16  
044324839

## Wiring diagram / Schaltplan

See the spare parts list for individual fuse part nos.

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**General information:**

The low-voltage fuse board protects all currents from the power transformer to the power output boards.

**Never use any fuse values or types other than those indicated**

The following table shows a list of fuses, each fuse being allocated to the part it protects.

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Keyboard - minimally equipped /  
Tastatur, min.


Abb. 3-5-1

Keyboard - fully equipped /  
Tastatur, max.


Abb. 3-5-2

**General information:**

The power output board controls the following devices:

- Chamber illumination
- Heating - Screen
- Heating - Window
- Heating - Tray
- Heating - Pipe defrost
- Solenoid valve MV1 (German abbr. MV1) Chamber
- Solenoid valve MV2 (German abbr. MV2) Shut-off valve - specimen
- Solenoid valve MV3 (German abbr. MV3) Control valve - specimen
- Fan control
- Heating cartridge rectifier
- Solid state relay compressor (German abbr. 'HBR')

Only the power output board for the 100V mains power input version has a different circuit. The difference lies in the control electronics for the solid state relay, which controls the compressors.

Following please find a description of the fuses for the power output board and for the protected devices.

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



Abb. 3-7-1

**Solder straps of drive boards and CM3050 S (0443 34561)**

These drive electronics have a total of 10 solder straps (denomination: SB1 ... SB10).  
For an overview, please see the tables which follow.  
Standard configurations (actual state of manufacturing) are underlined.

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| Terminal | Name of signals |             | Functions                          | Electrical characteristics      |
|----------|-----------------|-------------|------------------------------------|---------------------------------|
| CN1      | +24V            | Power input | ●Connect + and GND wires of 24V DC | 24V ±10%<br>CSK54: 1.3A maximum |
|          | GND             |             |                                    |                                 |

## Adjusting the driver output current CSK54

Procedures for lowering the output current to reduce the heat generation and vibration.

## Connections


Abb. 3-7-1-4

## Adjusting the running current CSK54

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
Abb. 3-7-1-5

## Adjusting the stop current CSK54

## Justierung des Ruhestroms CSK54


Abb. 3-7-1-3

## Precautions on installation

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Install from dust, oil mist, salt or corrosive gas???

Install from dust, oil mist, salt or corrosive gas???


Abb. 3-7-1-6

## Trouble shooting

When the stepping motor is not functioning properly, perform the following checks and the following results.

| Problem                                                              | Check Points                                                                                                | Measures                                                                                                                                                                           |
|----------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Motor is not energized.<br>(The motor shaft rotates easily by hand.) | ● Power Supply                                                                                              | ● Check that the power is connected                                                                                                                                                |
|                                                                      | ● All Windings Off input                                                                                    | ● When the All Windings Off input signal is "L" level (photocoupler On), the motor cease to be energized (has no holding torque)                                                   |
|                                                                      | ● motor                                                                                                     | ● Check that the motor and driver are connected properly                                                                                                                           |
|                                                                      | ● RUN and STOP potentiometers                                                                               | ● These potentiometers are used to adjust the output current to motor<br>(If they have been turned too far down, return them to their factory settings and then check the results) |
| Motor does not rotate<br>Motor rotates in the opposite direction     | ● Check the above points                                                                                    |                                                                                                                                                                                    |
|                                                                      | ● Pulse input or direction Input                                                                            | ● Check the connections, voltage and waveform of pulse signals                                                                                                                     |
| Motor is not functioning properly.                                   | ● Check the above points                                                                                    |                                                                                                                                                                                    |
|                                                                      | ● Are the motor and the load properly centered?<br>Is the load too large?                                   | ● Re-tighten the coupling screws or check with the load disengaged                                                                                                                 |
| Motor does not move far enough                                       | ● Does the actual motor step angle conform with the motor step angle required by the device?                | ● Check the setting of the step angle switch on the driver                                                                                                                         |
|                                                                      | ● Are the pulse generator settings for the input pulse number appropriate for the amount of motor movement? | ● Check the settings                                                                                                                                                               |

| Problem                                                                                                                  | Check Points                                                               | Measures                                                                                                                                                            |
|--------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The motor loses synchronism during acceleration (or during operation)                                                    | ● Is the All Windings Off signal input ?                                   | ● Disable the input                                                                                                                                                 |
|                                                                                                                          | ● Is the starting pulse too high ?                                         | ● Lower the rate and check the results                                                                                                                              |
|                                                                                                                          | ● Is the acceleration / deceleration time too short ?                      | ● Lengthen the time and check the results                                                                                                                           |
|                                                                                                                          | ● Is there enough current capacity in the power supply ?                   | ● Change to a power supply with a capacity of 2.1A or more and recheck (Then select the right power supply capacity)                                                |
|                                                                                                                          | ● Is there any effect from external electrical noise ?                     | ● Check the motor movement independently, without operating any other apparatus which could be potential sources of noise                                           |
| There is excessive vibration                                                                                             | ● There may be excessive motor output torque                               | ● Try reducing the motor running current                                                                                                                            |
|                                                                                                                          | ● Try changing the pulse rate                                              | ● If the vibration is reduced after changing the pulse rate , the problem might lie in the resonance of the motor Try changing the pulse rate or step angle         |
| The motor or the driver is excessively hot. (The temperature should be less than 100°C at motor case and 85°C at driver) | ● The motor has been operating for too long                                | ● Shorten the motor operating time or lengthen its rest time                                                                                                        |
|                                                                                                                          | ● Stop potentiometer setting is high                                       | ● Lower the current setting at motor current cut back                                                                                                               |
| The automatic current cutback function does not work                                                                     | ● Is the stop potentiometer in the maximum position ?                      | ● Current cannot be lowered when this potentiometer is in MAX position. Turn this potentiometer to the left , (Adjust to the optimal value by making reference page |
|                                                                                                                          | ● After conclusion of the pulse, is the pulse signal returned to H level ? | ● When the pulse signal is maintained at the "L" level (photocoupler is on) the current cannot be lowered. Be sure to return the pulse signal to the "H" position.  |

## Specifications

|             |                        |                                                   |                                                                                                                                                                                                    |
|-------------|------------------------|---------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Unit Model  |                        | Single Shaft<br>Double Shaft                      | CSK 545 - NATE<br>CSK 545 - NBTE                                                                                                                                                                   |
| Motor Unit  | Motor Model            | Single Shaft<br>Double Shaft                      | PK 545 - NA<br>PK 545 - NB                                                                                                                                                                         |
|             | Maximum Holding Torque | (N x m)<br>(oz - in)                              | 2,4<br>33                                                                                                                                                                                          |
|             | Rotor Inertia          | (kg x m <sup>2</sup> )<br>(oz - in <sup>2</sup> ) | 68 x 10 <sup>-7</sup><br>0,37                                                                                                                                                                      |
|             | Resistance Per Phase   | (Ohm/phase)                                       | 2,8                                                                                                                                                                                                |
|             | Mass                   | (kg)<br>(oz)                                      | 0,35<br>12,4                                                                                                                                                                                       |
|             | Step Angle             |                                                   | 0,72                                                                                                                                                                                               |
|             | Insulation Resistance  |                                                   | 100MOhm or more under normal ambient temperature and humidity when the megger reading between the windings and the frame is DC500V.                                                                |
|             | Dielectric Strength    |                                                   | Under normal ambient temperature and humidity, sufficient to withstand 1kV at 50Hz applied between the windings and the frame for one minute following a period of continuous operation.           |
|             | Insulation Class       |                                                   | Class B (130°C)                                                                                                                                                                                    |
| Driver Unit | Driver Model           |                                                   | CSD5807N - T                                                                                                                                                                                       |
|             | Voltage                |                                                   | 24V DC +/- 10%<br>1,3A maximum                                                                                                                                                                     |
|             | Output Current         |                                                   | 0,1 ~ 0,75A / phase                                                                                                                                                                                |
|             | Excitation Mode        |                                                   | Full step: 0,72° / step<br>Half step: 0,36° / step                                                                                                                                                 |
|             | Input / Output Signals |                                                   |                                                                                                                                                                                                    |
|             | Pulse Input            |                                                   | 1 - pulse input<br>Photocoupler input<br>Input Impedance 220Ohm, Input current 20mA maximum<br>Pulse width 5μ sec. minimum, Pulse rise / pulse fall time 2μ sec. maximum<br>L: 0 ~ 0,5V, H: 4 ~ 5V |
|             | All Windings Off input |                                                   | Photocoupler input<br>Input Impedance 220Ohm, Input current 20mA maximum<br>L: 0 ~ 0,5V, H: 4 ~ 5V                                                                                                 |
|             | Full / Half Step input |                                                   |                                                                                                                                                                                                    |

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
Abb. 3-7-1-7




Abb. 3-7-1-8

## Feed unit / Zustelleinheit

Make sure to  
install in correct  
position!  
Einbaulage  
beachten!

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









Circuit diagram, Schaltplan

U2-11

### Adjustment Sectioning Motor Slow / Fast

Tool: Screw driver

#### I. Zero drift

Set the sectioning motor speed to rapid prior to initialising the instrument.

Select continuous stroke mode of operation; set a large sectioning window and start the instrument.

Set the speed potentiometer to zero and check at the sectioning motor whether all movement has actually completely ceased.

In case the motor is still moving, bring it to a halt by adjusting potentiometer R39.

#### II. Halt point

Set the speed potentiometer between zero and the first index line and adjust R39 so the motor rotates at a very slow speed.

Set the speed potentiometer to zero and, at the potentiometer, check the halt point.

Perform an instrument RESET, or rather switch the instrument off and back on.

#### III. Adjustment sectioning motor-slow

Prior to initialising the instrument, select 'Sectioning motor – slow'.

Set zero drift and halt point the same way as described under I. and II.

The adjustments have to be done with potentiometer R38.

Once the adjustment procedure is finished, safeguard the potentiometer with sealing lacquer.

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
Abb. 3-10-1

**Note**

Each menu point can be switched ON / OFF and its value can be changed by pressing the "DEFROST KEY".

When scrolling through the menu, some consumers are automatically switched off, if they are called up for the second time.

(Exception: Specimen ... menu points are switched off immediately once scrolling is continued).

To finish the service routine, press "SCROLL".

Once the service routine has been finished, the instrument initializes itself (reset).

**Note on the Service Routine:**

The service routine is a tool which enables an expert to perform a simple check on the instrument without any other tools being required. Together with the corresponding control LEDs situated on each board, the service routine facilitates an easy-to-do service.

**Attention!**

The service routine is not a diagnostic routine which is executed automatically. To run the service routine, the user needs to have a corresponding level of technical knowledge. Please note that operation errors in running the service routine can even cause damage to the instrument.

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Press key "KEY", "ARROW UP" and "LIGHT"

(does not work in standby mode)

The following menus can be chosen by pressing the "ARROW" keys.

|                                                                                   |                                                      |
|-----------------------------------------------------------------------------------|------------------------------------------------------|
| 1. POT = XXX<br>BRAKE = ZZ                                                        | XXX = tandempot is displayed<br>ZZ = set brake value |
| 2. WDOG                                                                           | ON/OFF                                               |
| 3. Specimen fast FWD                                                              | ON/OFF                                               |
| 4. Specimen fast BWD                                                              | ON/OFF                                               |
| 5. Specimen slow FWD                                                              | ON/OFF                                               |
| 6. Specimen slow BWD                                                              | ON/OFF                                               |
| 7. Cutmotor (additional RUN ENABLE-key necessary)                                 | ON/OFF                                               |
| 8. set default clock-default adjustment                                           | ON/OFF                                               |
| 9. Valve (Cabin) (defrost valve)                                                  | ON/OFF                                               |
| 10. Magn.Valve Odj. (object cooling)                                              | ON/OFF                                               |
| 11. Reg.Valve Obj.                                                                | ON/OFF                                               |
| 12. Runtime... (in hours) delete with ON<br>reactive with OFF, resp. by scrolling | ON/OFF                                               |
| 13. Heating (Object)                                                              | ON/OFF                                               |
| 14. Heating (window, pane)                                                        | ON/OFF                                               |
| 15. Heating (Cabin)                                                               | ON/OFF                                               |
| 16. Compressor (chamber cooling)                                                  | ON/OFF                                               |
| 17. Compressor (object cooling)                                                   | ON/OFF                                               |
| 18. Light                                                                         | ON/OFF                                               |
| 19. Display Light                                                                 | ON/OFF                                               |
| 20. Button LEDs                                                                   | ON/OFF                                               |
| 21. Trim LEDs                                                                     | ON/OFF                                               |
| 22. Clutch                                                                        | ON/OFF                                               |
| 23. Object Position XX                                                            | 0 – 255 Typical value: 215 + 3                       |
| 24. Cab. Temp.                                                                    |                                                      |
| 25. Temp.                                                                         |                                                      |
| 26. ADC Cab. Ref.                                                                 |                                                      |
| 27. ADC Obj. Ref                                                                  |                                                      |

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
**Important note!**

The operating instructions below are valid only for the instrument Leica CM 3050 S. They supplement chapter 7.5 of the standard Leica CM3050 instruction manual, explaining the new / additional functions of the Leica CM3050 S unit. All references to figures, page numbers and / or headings relate to the instruction manual version 1.3-07/96.

**1. Mains switch**

The instrument is no longer switched on/off via the mains switch shown in Fig. 21.1 (which has been eliminated in the new Leica CM3050 S) but via switch (1), which combines the function of ON/OFF switch and automatic cutout for the instrument. eingeschaltet. The new layout looks as shown on the left.

**2. Control panel 1**

Control panel 1 as shown in Fig. 22.1 will not be modified as far as shape and key functions are concerned. New is the instrument name **LEICA CM3050 S**, as well as a number of new display indications.

**3. Control panel 2**

Module B of control panel 2 shown in Fig. 24.1 has been changed.

Modules A and C remain the same, i.e. are identical with those in the Leica CM 3050 unit.

Module B has been modified as follows:


**1. Trimming**

The  button has been replaced by the  button.


To switch on press , TRIM (2) illuminates. The illuminated bar (3) indicates the selected trimming thickness in  $\mu\text{m}$ .


**2. Selection keys**

The arrow buttons used in the Leica CM3050 have been replaced by the  and  buttons. Pressing  decreases the selected value (see display indication), pressing  increases the selected value.

## 3. Control panel 1 Display (Fig. 23.1, 23.2)

'Preset counter' deactivated

**Upper line:**

Immediately after switching on, the instrument name CM3050 S and the installed software version are displayed for a short time.

**Lower line:**

„MOT“ for „motorized sectioning on“ is not displayed any longer. The 8 digits on the right of the lower line (see marked area) are now used to display a number of additional parameters. All other settings shown in Fig. 23.1 and 23.2 remain the same.

By pressing  a number of parameters are displayed one after another in the lower line (sequence as shown on the left). By pressing  the same sequence of parameters is shown, however in reverse order.

**1. Section thickness**

Display of selected section thickness in  $\mu\text{m}$ . Select section thickness with  /  buttons of control panel 2 while sectioning mode is activated (Indication TRIM --> OFF).

**2. Section counter**

Counts the sections carried out. Press  to set back to „0“.


**3. Section thickness totalizer**

Displays the total thickness of all sections carried out. Indication in  $\mu\text{m}$ . Press  to set back to „0“.


**4. Preset Counter**

To carry out a preselected number of sections. When sectioning is started, the preset counter starts counting backwards. When dn=000 is reached, sectioning in continuous mode stops. To set back to the preselected value, press button  and  simultaneously. Press  to set back to „0“.

The preset counter can be deactivated. To activate, display 'preset counter' pressing the menu button (see point 4) then press  to activate. Prior to starting the preset counter, a value (= pre-selection of number of sections) must be entered.

**4. Control panel 1****Menu button**

4


The menu symbol itself (4) remains the same. However, the sequence of indication of the individual values has been changed.

By pressing the menu button, the different parameters are displayed in the following order (press menu button once per parameter):

1. Chamber temperature
2. Specimen temperature
3. Preselected value of preset counter  
(Press  to activate / deactivate preset counter)
4. Maximum (= lowest) temperature of specimen head
5. Time (real time)
6. Defrost time (= time when defrost cycle starts)
7. Duration of defrost cycle (Chamber)

**5. Control panel 2**

 /  - Tasten

When the instrument is in sectioning mode, (indication TRIM --> OFF!) display indication will go to sectioning thickness selection by pressing  or .

Displayed is always the section thickness value which has been selected last.

Press  or  to select a new section thickness value.

Selectable section thickness selection range:

- |     |   |        |           |        |
|-----|---|--------|-----------|--------|
| 0   | - | 2 µm   | in 0.5-µm | steps, |
| 2   | - | 10 µm  | in 1-µm   | steps, |
| 10  | - | 20 µm  | in 2-µm   | steps, |
| 20  | - | 60 µm  | in 5-µm   | steps, |
| 60  | - | 100 µm | in 10-µm  | steps, |
| 100 | - | 300 µm | in 50-µm  | steps. |

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


4.5 (Box)



Abb. 960919k1

## Disassembly / assembly of a compressor (cryochamber refrigeration)

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**Introduction**

Before starting to perform repair work on a refrigeration system, make sure that the premises where the repair will take place are adequately ventilated. An essential part of any good assembly / repair work is adequate preparation.

We also take it for granted, that soldering is done with protective gas (flooding with the pipe with nitrogen or other protective gas type - no pressure build-up).

1. Install service valve at both pressure and suction line.
2. Recover refrigerant from refrigerating system and perform pressure compensation.
3. Exchange the compressor.
4. Exchange the filter-drier.  
**Attention!** In all refrigerating plants filter-driers are used in the liquid pipe. Such as the compressor can be regarded as the 'heart' of the refrigerating system, the filter-drier can be considered the 'liver' of the refrigerating system. An arrow on the filter-drier indicates the direction of flow. When installing a filter-drier make sure to not open it until just prior to installation into the ready-for-service refrigerating system. That way you ensure that the filter-drier does not absorb any ambient humidity prior to installation. The filter-drier absorbs the following substances from the refrigerant circulating refrigeration cycle:

- a) solid impurities (e.g. chips, cinder, oleoresin, tinder)
- a) Water
- c) Acid

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

Abb. 960919k1

5. Evacuate the refrigeration cycle: suction and pressure pipes at least 3.0 hours and using a two-step evacuating pump.  
Once the evacuation procedure is finished, there should be a vacuum of at least 0.3 mbar.

**Attention!**

**If the vacuum does not remain stable after evacuation, possibly there has to be done a pressure check. Only dry nitrogen should be used to do the pressure check. Make sure to attach a reducing valve to the nitrogen cylinder valve prior to filling nitrogen into the refrigerating system. Search for leaks very carefully.**

6. Add refrigerant (80% of the quantity indicated on the instrument name plate).

**Attention!**

**Only use liquid refrigerant to fill the refrigerating system.**

7. Make sure to fill the liquid refrigerant directly into the liquid pipe.
8. Switch the compressor on.
9. Adjust the refrigerant cycle. Overheating of the gilled evaporator has to be 2 - 3 K. Adjust overheating by adding or removing refrigerant.
10. Seal off the service pipe, remove the service valves and solder up the opening in the copper pipe.

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**Trouble-shooting refrigerating systems**

High-quality repairs must be preceded by a systematic trouble-shooting procedure!

The most frequent malfunctions can be roughly categorized as follows:

- electrical malfunctions, including safety shut-downs.
- mechanical defects of compressor and functional elements
- insufficient refrigeration

The table below is therefore only a short set of instructions on how to proceed.

**Malfunction 1: Compressor does not start**

| Possible cause                                                              | How to verify the cause                                               | Remedy                                                 |
|-----------------------------------------------------------------------------|-----------------------------------------------------------------------|--------------------------------------------------------|
| Mains switch off                                                            | Check switch position                                                 | Select correct switch position (ON)                    |
| Fuse activated / blown                                                      | Check fuse (visual inspection or with load testing device)            | Replace the fuse or push it back in                    |
| Power supply to control box interrupted                                     | Use load testing device to verify voltage feed at input clamps        | Have damage in mains or lead removed by electrician.   |
| Mechanical start strain relief missing or defective                         | Visual inspection, functional check of solenoid valve and check valve | Upgrade instrument or exchange defective part          |
| Motor coil defective                                                        | Check resistance of coils                                             | Exchange motor or entire compressor                    |
| Mechanical compressor damage such as broken connecting rod, damaged bearing | Visual inspection                                                     | Exchange compressor                                    |
| Temperature controller defective or switched off                            | Electrical check                                                      | Exchange compressor or temperature controller (Klixon) |

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**Malfunction 2. Compressor starts, but stops again after a short time**

| Possible cause                                                        | Verifying possible causes                                                                                         | Remedy                                                                                                                 |
|-----------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------|
| Condenser pressure too high<br>(Suction pressure normal or increased) | Measure condenser pressure and ambient temperature, check condenser, check dirt accumulation on condenser surface | Improve ventilation of installation site, clean condenser, evacuate refrigeration cycle, exchange defective components |
| Condenser pressure too high<br>(Suction pressure too low)             | System clogged, e.g., filter-drier, capillary tube<br>Expansions valve                                            | Replace defective parts                                                                                                |

**Malfunction 3. Compressor running, but refrigeration is insufficient**

| Possible cause                                                                                                 | Verifying possible causes                                                                         | Remedy                                                                                  |
|----------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|
| Condenser wear and tear or mechanical damage. Valve reed or cylinder head gasket defective. Piston rings worn. | Check cylinder head temperatures. Check pressure ratio.                                           | Remove compressor.                                                                      |
| Insufficient quantity of refrigerant.                                                                          | Measure suction pressure<br>Check frost build-up on evaporator                                    | Find leaks and eliminate.<br>Add refrigerant.                                           |
| Secondary lack of refrigerant<br>Filter-drier clogged.                                                         | Suction pressure too low.<br>Insufficient frost build-up on evaporator.<br>Excessive overheating. | Remove blockage.<br>Replace defective components.                                       |
| Fan defective.                                                                                                 | Measure suction pressure.                                                                         | Replace fan motor.                                                                      |
| Evaporator heavily iced .                                                                                      | Check if defrosting functions properly, check location of instrument                              | Defrost evaporator.<br>Eliminate defect(s) in automatic defrost system, change location |

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
Abb. scan3a

**Assembly / disassembly of the control valve**

1. If there is any visible damage on the sealing surface of the housing (1) and/or on the sealing surface of the solenoid valve housing (4), we recommend to exchange the entire control valve.
2. If the sealing surfaces of the solenoid valve housing are not defective, it is sufficient to exchange only the screwed valves.

**Attention! The screw valve of control valve (2) has to be tightened at a torque of 20 Ncm.**

3. If the entire control valve unit or the screwed valves have to be exchanged, it is necessary to open the refrigerating system (i.e. the refrigerant must be recovered).
4. The control valve unit assembly is soldered into both openings of the output side of the specimen head (see chapter 4.4).
5. Adjust the refrigerating system (see chapter 4.4).

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
Abb. 311296k1

**Installing the specimen head**

1. Inside the cryochamber, insert the specimen head from above into the corresponding opening and secure with four screws (2).
2. Fit copper pipes (see chapter 4.2) and copper pipe output of the specimen head into each other and solder up.
3. Fit capillary tube (1) into the copper pipe of the control valve and solder it in.
4. Fit the copper pipe input of the specimen head into the copper pipe of the control valve and solder it in.
5. Install a service valve on both service pipes (pressure and suction pipe).
6. Evacuate refrigerating system for at least three hours with a two-stage vacuum pump.
7. Fill refrigerant into refrigerating system, see adjustment instructions.
8. Press shut both service pipes right before the service valves and then remove the service valves. Solder up the openings.

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
Abb. 311296k1

**Instructions for filling and adjusting the specimen cooling**

1. Turn the instrument on (mains switch).
2. Turn off both compressors in standard operation mode:  
Press the MENU button until SET TEMP OT—°C and SET TEMP CT—°C then appear. Deactivate the specimen cooling system as well as the chamber cooling system pressing the KEY button (compressors off).  
Two hyphens are indicated instead of asterisks.
3. Turn off instrument via mains switch.
4. Recover refrigerant, if necessary, perform repair, replace filter-drier.
5. Turn on the instrument mains switch.
6. Activate the service routine by pressing the KEY, ARROW UP and LAMP buttons simultaneously.
7. Press the ARROW buttons until REG. VALVE OBJ. is displayed.
8. Press the MANUAL DEFROST BUTTON to activate REG. VALVE OBJ. (= valve open).
9. Evacuate refrigeration system for at least 3.0 hours.
10. Add Frigen (refrigerant) (80% of the quantity indicated on the instrument nameplate).
11. Press the MENU button to leave the service routine. The instrument is then automatically initialized (reset).
12. Activate both compressors via the standard operating mode (KEY button - see point 2.)
13. Set chamber temperature to -35 °C and specimen temperature to -50 °C and wait until the instrument has cooled down to the selected chamber temperature.
14. Check the lower limit temperature of the specimen cooling system. For that purpose the chamber temperature needs to be stable. Depending on the ambient temperature, the temperature of the specimen head varies as indicated in the enclosed diagram.  
Pull off the solenoid of the control valve to check these temperature values.  
Tolerance: +/- 2 K  
Example: At an ambient temperature of +26 °C, the specimen head temperature should be -50°C +/- 2K.
15. If the specimen temperature is somewhere in the range of -53 °C ... -56 °C, there is an insufficient quantity of refrigerant: add refrigerant in 5-g steps.

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
Abb. 311296k1

If the specimen temperature is somewhere between -47 °C.... -44 °C, too much refrigerant has been added - remove refrigerant in 5-g steps.

**Important!**

**Refrigerant R404A (Frigen) MUST be in a liquid state when being added into the refrigeration cycle.**

**After each step (i.e. adding or removing refrigerant), a waiting period of approx. 5 minutes must be observed.**

16. Reinstall the solenoid of the control valve.

Set specimen temperature to -10°C.  
Wait approx. 5 min; then check the actual temperature. A tolerance of +/- 2 K is acceptable.

**Note:**

**If the actual temperature is lower than the set value (i.e. colder) —> Add refrigerant.**

**If the actual temperature is higher than the set value (i.e. warmer) —> Remove refrigerant.**

**Important!**

**Refrigerant R404A (Frigen) MUST be in a liquid state when filling.**

**After each step (i.e. adding or removing refrigerant), a waiting period of approx. 5 minutes must be observed.**

17. Once the correct temperature value has been obtained, proceed to check the cooling-down time.

Cooling-down time:

|                    |       |
|--------------------|-------|
| From OT to OT Set  | max.  |
| -10 °C      -45 °C | 140 s |

Tolerance:

Temperature: +/- 1K

Time: +10 s

18. If the period of time measured is shorter than required and the temperature is e.g. -47 °C...-43 °C and higher (warmer)  
-> remove refrigerant in steps of 5g  
If the period of time measured is longer than required  
-> refill refrigerant in in steps of 5g

**Important!**

**Refrigerant R404A MUST be in a liquid state when filling.**

**It is essential to wait approx. 5 min after every modification, i.e. after extracting or refilling refrigerant.**

19. Press shut service pipes, remove service valve, solder up openings.

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
Abb. 4-4-1

Legend:

- 1 Compressor
- 2 Condenser
- 3 Fan motor
- 4 filter-drier
- 5 Gilled evaporator capillary
- 6 Gilled evaporator
- 7 Wall evaporator capillary
- 8 Wall evaporator and quick-freeze shelf
- 9 2/2 Solenoid valve with magnet coil
- 10 Liquid collector
- 11 Filling and recovering pipe
- 12 Check valve
- 13 Insulation material


Abb. 960710k1

| Temperature sensor pos | Measurement point        | should be                                                        | actual | Comments                                                                                                                                          |
|------------------------|--------------------------|------------------------------------------------------------------|--------|---------------------------------------------------------------------------------------------------------------------------------------------------|
|                        | Room temperature         | + 23°C                                                           |        |                                                                                                                                                   |
|                        | Box temperature          | - 25°C                                                           |        |                                                                                                                                                   |
| 3                      | tv 1                     | 0°C ... + 5°C                                                    |        | measurement pos. approx. 20cm in front of compressor (turn of the unit, 15min. measurement to be done within 3min. after starting the compressor) |
| 4                      | tv 2                     | + 85°C ... + 95°C                                                |        | measurement pos. approx. 10cm behind compressor                                                                                                   |
| 5                      | tc 2                     | + 28°C ... + 35°C                                                |        | measurement on bottom tube coming from condensor                                                                                                  |
| 1                      | to 1 (gilled evaporator) | - 42°C ... - 47°C                                                |        | measure on evaporator in chamber, pipe on left side above drip pan rear side                                                                      |
| 2                      | to 2 (gilled evaporator) | - 38°C ... - 44°C                                                |        | measure on evaporator in chamber, pipe on left side above drip pan front side                                                                     |
|                        | tb 1                     | + 62°C ... + 68°C                                                |        | measure on tube going into defrost valve in front of solder point                                                                                 |
|                        | to (freezing stage)      | - 38,5°C ... - 44°C                                              |        | measurement position on the freezing stage directly                                                                                               |
|                        | to 1 (gilled evaporator) | + 55°C ... - 44°C                                                |        | during defrost period 9min.                                                                                                                       |
|                        | po                       | 1.1 ... 1.0bar abs<br>to = - 46°C ... -48°C<br>temp. scale R404A |        | hook up gauges on service valve low side                                                                                                          |
|                        | Pressure balance         | 7.3bar                                                           |        | if the instrument is not in function                                                                                                              |
| 6                      | Compressor head          | + 53°C ... + 60°C                                                |        |                                                                                                                                                   |

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Abb. Tab1

Legend:

- 1 Compressor
- 2 Condenser - specimen
- 3 Fan motor
- 4 Filter-drier
- 5 2/2 Solenoid valve
- 6 2/2 Solenoid valve with magnet coil
- 7 Capillary tube
- 8 Flexible metal tube
- 9 Specimen head
- 10 Heating cartridge
- 11 Heat exchanger
- 12 Filling and recovering pipe
- 13 Insulation material


Abb. 960710k2

| Temperature sensor pos | Measurement point    | should be                                                                           | actual | Comments                                                      |
|------------------------|----------------------|-------------------------------------------------------------------------------------|--------|---------------------------------------------------------------|
|                        | Room temperature     | + 23°C                                                                              |        |                                                               |
|                        | Specimen temperature | - 50°C                                                                              |        |                                                               |
|                        | tv 1                 | + 15°C ... + 25°C                                                                   |        |                                                               |
|                        | to const.            | - 50°C                                                                              |        | measured at specimen holder with temperature measuring device |
|                        | to const.            | - 10°C ... - 4°C<br>140s.<br>- 40°C ... - 5°C<br>100s.<br>- 50°C ... - 10°C<br>60s. |        | Display                                                       |
|                        | po                   | 0.7bar abs<br>to - 54°C<br>(temp. scale R404A)                                      |        | hook up gauges on service valve low side                      |
|                        | Pressure balance     | 7.0bar                                                                              |        | if the instrument is not in function                          |
| 8                      | Compressor head      | + 50°C - + 55°C                                                                     |        |                                                               |
| 7                      | tc 2                 | + 27°C - + 30°C                                                                     |        |                                                               |

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**Spare parts can be obtained from your local Leica office or from**

Leica Microsystems Nussloch GmbH  
Postfach 1120  
Heidelberger Str. 17-19  
D-69222 Nussloch  
Germany  
Telephone: 49 6224 1431 163  
Facsimile: 49 6224 1431 253

**Technical Service - Hotline**

There are two service hotlines for customer service and requests:

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Telephone: + 6224 143 229 for electronical requests  
Telephone: + 6224 143 219 for refrigerational and mechanical requests  
Facsimile: + 6224 143 199 for electrical and mechanical requests

The agencies are obliged to send complete failure statistics (monthly) to Nußloch for the first year after instruments have been launched onto the market. The results and necessary amendments can be included immediately in the production.

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


| No. | Part Number | Qty. | Description                       | Remarks |
|-----|-------------|------|-----------------------------------|---------|
| 1   | 0443 33466  | 1    | Keyboard, assy.<br>Keyboard, kpl. |         |

| No. | Part Number | Qty. | Description                                                         | Remarks |
|-----|-------------|------|---------------------------------------------------------------------|---------|
| 1   | 0443 33507  | 1    | Keyboard CM3050 S, fully equipped<br>Flächentastatur CM3050 S max.  |         |
|     | 0443 33506  | 1    | Keyboard CM3050 S, minimally equipped<br>Flächentastatur 3050S min. |         |

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